

HEINE QUALITY
MADE IN GERMANY



LED HO

LED NOW IN HEINE QUALITY.

Quality in a new light
: **HEINE OMEGA[®] 500**
Indirect Ophthalmoscope





Our Best Indirect Ophthalmoscope
is now twice as good:
With LED- and XHL-Technology

The focus remains providing the ultimate in Quality, Performance, and Flexibility.

Ultimately you are investing in an indirect ophthalmoscope. For HEINE, the quality of your indirect system begins with the optics. The OMEGA 500 offers adjustments and flexibility not available in any other instrument, and utilises the most durable materials to provide you with the performance and reliability you demand.



STEREO VIEWS REGARDLESS OF PUPIL SIZE

The OMEGA 500 gives you total freedom to adjust your optical alignment to maximise your view in all pupils. The exclusive synchronised convergence and parallax adjustment system is free of presets, allowing stereo views in pupils ranging from 1 mm through 10 mm.



PERFORMANCE AND FLEXIBILITY

With the OMEGA 500, all you need is right at your fingertips. Three spot sizes and a diffuser are integrated, along with three filters: yellow filter for patients who are sensitive to light, red-free filter, and a cobalt blue filter. The tiltable mirror allows you to adjust the light beam + or - 4 degrees.



MAINTENANCE FREE

The 100 % guaranteed dustproof construction and mounting of the optics on an aluminum frame make the OMEGA 500 maintenance free and durable. It is so reliable that we have extended its guarantee to 5 years.



DESIGNED FOR COMFORT AND ERGONOMICS

The unique headband of the OMEGA 500 is perfectly balanced and comfortable due to its strong and malleable materials and exclusive adjustment possibilities. Height, distance, and angle positioning for your optics, combined with independent interpupillary adjustment and ergonomic control placement allow for a customised fit and intuitive operation.





XHL or LED?

STILL THE STANDARD:
HEINE XHL-TECHNOLOGY

XHL

RELIABLE, ACCURATE, FAMILIAR: HEINE XHL-TECHNOLOGY

HEINE XHL Xenon-Halogen-Technology has set the standard for the evaluation of new light technologies because it offers the perfect colour rendering that health professionals have examined the eye with and depended on for decades.

But now there is another illumination option that meets this standard: our new LED-Technology. Now you have a choice with no compromises: tried and true XHL, or LED in HEINE Quality.



FUTURE CHANGE OF MIND NOT EXCLUDED

Times are changing. Opinions too. And sometimes it's just the desire for new things that prevails. It doesn't matter if it's a brand new or an older OMEGA 500; you have the ability to easily change between Halogen- and LED-illumination. For more information on upgrading from XHL, contact HEINE or your authorised distributor.

Which is right for you?

SEE THE EYES IN A NEW LIGHT: WITH HEINE LED



EXACT BRIGHTNESS CONTROL

For proper brightness control and flexibility we've adapted the rheostat to accommodate for the unique characteristics of LED- and XHL-Technologies. The HC50 L can be used in conjunction with the LED- and the XHL-system without any alterations or switches.



HOMOGENEITY

To be able to discern even the slightest alterations of the retina, a homogeneous, distortion-free illumination of the field of view is important. The optics and illumination system of the OMEGA 500 have always been known for this feature, and LED^{HQ} is no exception.

**CRI
90**

COLOUR RENDERING

The Colour Rendering Index (CRI) indicates how exactly colours are reproduced in comparison with daylight. The closer this value is to the maximum of 100, the more accurate the lighting system is. The OMEGA 500 with LED produces an excellent CRI value of 90, and is strong in its rendering of the colour red in particular, with a value of 75. Combined with a colour temperature of 4000K, this means that should you choose LED or opt to change from XHL- to LED-Technology, the colour rendering remains almost unaltered.



THERMAL MANAGEMENT

Proper thermal management is required for a LED illumination to ensure consistent light output and maximum lifetime. Our LED module is constructed of heat conducting materials to assure a fast dissipation of the heat. This guarantees a consistent brightness of the LED during the whole working life.

**LONG
LIFE**

LIFESPAN AND ON TIME

By utilising the latest LEDs and highest quality materials, our LED-Technology ensures a virtually unlimited working life (typ. 50,000 hours). Additionally, rechargeable batteries last longer due to the increased operating time and the reduced number of charging cycles.



LED NOW IN HEINE QUALITY.

GUARANTEE



HEINE makes no compromises in manufacturing high quality medical instruments. Our commitment to vertical integration in manufacturing means we control all aspects of our instrument quality, from utilising carefully selected, well-matched materials to a high level of hand assembly. This ensures that each HEINE instrument meets or exceeds all requirements in any medical environment. We stand behind this commitment with a 5 year guarantee!



No Detail is Overlooked

COMFORT BEGINS WITH THE FIT AND ERGONOMICS

WORKING IN PERFECT BALANCE

The aim of the OMEGA 500 headband: maximum wearing comfort. The inclined overband reduces downward force which relieves strain on your neck muscles. The strong, malleable material and padding form to your anatomy for a soft, secure fit without pressure points.

INTUITIVE OPERATION AND FEATURES

The ergonomically arranged controls for the selection of apertures and filters are equipped with soft-touch surfaces and are easily adjusted while wearing the instrument. Exclusive locking mechanisms will keep your controls in the desired position, protecting against accidental changes, but will disengage without taking off the instrument. The cables are integrated in the headband and contribute to a safe handling and operation.

FLIP-UP SYSTEM FOR EASY COMMUNICATION

The exclusive locking flip-up system of the headband has been developed to allow you to securely raise and lower your optics. The optics can be raised and securely fixed without changing the settings, for an easier communication with the patient.

RHEOSTAT POSITIONING

Why reach across your body to adjust the brightness control of your indirect ophthalmoscope? With the OMEGA 500 you will never have to, because the HC50L Rheostat can be mounted on either the left or right side of the indirect ophthalmoscope. The change is quick and easy and can be done by the users on their own.

Power Your Way

FLEXIBLE AND MODULAR POWER SYSTEM

A SYSTEM FOR EVERY ENVIRONMENT

Because of the quality of the OMEGA500, it often accompanies its owner for a life time. From school and training, through to a practice or hospital environment. With the most complete and flexible power source offering, your OMEGA500 will adapt to any environment. If your needs or requirements change, the power sources are modular and completely interchangeable.

mPACK PORTABLE RECHARGEABLE BATTERY

Rechargeable battery that you can wear clipped on your belt or in your pocket.

- Up to 20 hours of ON time with LED illumination.
- Minimum 4 hours of ON time with XHL illumination.
- 2 hours quick charge time.
- Exclusive 5 step battery level display tells you the remaining battery charge as well as the charge status.
- Our exclusive mobile charger allows you to power your instrument and charge it at the same time.
- For office based environments we offer the wall mountable EN50 Transformer.

mPACK UNPLUGGED HEADBAND-MOUNTED RECHARGEABLE BATTERY

Complete freedom, no cables, and no compromises regarding performance.

- More than 8 hours of ON time with LED illumination.
- Minimum 2 hours of ON time with XHL illumination, depending on the brightness setting.
- 2 hours quick charge time.
- Exclusive 5 step battery level display tells you the remaining battery charge as well as the charge status.
- Our exclusive mobile charger allows you to power your instrument and charge it at the same time.
- For office based environments we offer the Wall Transformer EN50 UNPLUGGED with automatic ON/OFF function and the ability to charge a spare battery.

TABLE-TRANSFORMER EN50

If a mobile operation isn't necessary, the OMEGA500 can be operated with the table or wall mountable EN50 Transformer.

OMEGA500 UNPLUGGED Kit with Wall Transformer EN50 UNPLUGGED

Kit containing the OMEGA500 with LED- or XHL-Technology, the HC50 L Headband Rheostat, the Wall Transformer EN50 UNPLUGGED and 2 Headband-mounted Rechargeable Batteries mPack UNPLUGGED.

LED-Illumination: C-008.33.536

XHL-Illumination: C-004.33.542

OMEGA500 UNPLUGGED Kit with Headband-mounted Rechargeable Battery mPack UNPLUGGED

Kit containing the OMEGA500 with LED- or XHL illumination, the HC50 L Headband Rheostat, a Headband-mounted Rechargeable Battery mPack UNPLUGGED and the corresponding Plug-in Transformer.

LED-Illumination: C-008.33.535

XHL-Illumination: C-004.33.541

OMEGA500 Kit with mPack Rechargeable Battery

Kit containing the OMEGA500 with LED- or XHL-Technology, the HC50 L Headband Rheostat, the Rechargeable Battery mPack and the corresponding Plug-in Transformer.

LED-Illumination: C-008.33.533

XHL-Illumination: C-004.33.539





In 1946 our founder, Helmut A. Heine, started HEINE Optotechnik with one goal in mind; to build the highest quality diagnostic instruments available in the world. It was his belief that the highest quality instruments will always best support the health care professional in making the earliest, and most accurate possible diagnosis. For over 70 years, our customers have been able to rely on this promise and on their HEINE instruments for precise diagnostic performance, unique functionality and long lasting durability.

Today, as a global leader in primary diagnostic instruments with over 500 employees we remain committed to our founding principles. We are still a 100% family owned and managed company. We care deeply about our employees, our customers and the quality of our products. We continue to develop and manufacture HEINE instruments exclusively at our facilities in Germany, where we combine decades-long experience and craftsmanship with the most modern manufacturing technologies.

www.heine.com



We reserve the right to change specification without notice.
© 2018 HEINE Optotechnik GmbH & Co. KG. All rights reserved.

Your distributor

03/18. A-000.00.221 e

GERMANY

 HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7
82211 Herrsching
Tel. +49 (0) 81 52-38 0
Fax +49 (0) 81 52-3 82 02
E-Mail: info@heine.com
www.heine.com

NORTH AMERICA

HEINE USA LTD.
10 Innovation Way
Dover, NH 03820
Tel. (603) 7 42-71 03
Fax (603) 7 42-72 17
Toll Free (800) 367-4872
E-Mail: service@heine-na.com

AUSTRALIA

HEINE AUSTRALIA PTY. LTD.
Unit 9, 98 Old Pittwater Road
PO Box 7218 Warringah Mall
NSW 2100
Tel. +61 (0) 2-99 38 95 00
Fax +61 (0) 2-99 39 23 05
E-Mail: info@heine.com.au

SWITZERLAND

HEINE (Schweiz) AG
Tobeläckerstr. 9
CH-8212 Neuhausen am Rheinfall
Tel. +41 (0) 52-6 72 22 66
Fax +41 (0) 52-6 72 63 77
E-Mail: info@heineschweiz.ch

